



1 Cephaelin  
0.017 g/mL

2 Emetine  
0.017 g/mL

#### HPLC Conditions

Column CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm  
Mobile phase 2w/v% Sodium 1-Heptanesulfonate (pH4.0, CH<sub>3</sub>COOH) / CH<sub>3</sub>OH = 50 / 50  
Flow rate 0.5 mL/min  
Temperature 50 °C  
Detect UV 283 nm  
Inj. vol. 10 µL