



1	5'-CMP 100 µg/mL	2	5'-AMP 100 µg/mL	3	5'-UMP 100 µg/mL
4	5'-IMP 100 µg/mL	5	5'-GMP 100 µg/mL		

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d.×250 mm
Mobile phase	50mmol/L (NH ₄) ₂ HPO ₄ (pH3.0)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 260 nm
Inj. vol.	10 µL