



1	Uracil	2	Thymine	3	Guanine
	5 µg/mL		5 µg/mL		15 µg/mL
4	Cytosine	5	Adenine		
	30 µg/mL		15 µg/mL		

HPLC Conditions

Column	CAPCELL PAK SCX UG80 : 4.6 mm i.d.×250 mm
Mobile phase	200mmol/L NH ₄ H ₂ PO ₄ (pH3.50)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm
Inj. vol.	10 µL