



- | | | |
|--------------------------|---------------------------|--------------------------|
| 1 Prednisolone | 2 Flucinolone acetonide | 3 Hydrocortisone acetate |
| 4 Dexamethasone acetate | 5 Hydrocortisone butyrate | 6 Flucinonide |
| 7 Chlobetazol propionate | | |

HPLC Conditions

Column	CAPCELL PAK C18 AG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH ₃ OH / CH ₃ CN / H ₂ O = 40 / 23 / 10
Flow rate	1.0 mL/min
Detect	UV 243 nm