



1	Estriol	2	$\beta$ ]Estradiol	3	Estrone
	10 ppm		10 ppm		10 ppm
4	$\beta$ ]Estradiol-3-methylester				
	10 ppm				

#### HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH <sub>3</sub> CN / H <sub>2</sub> O = 55 / 45
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 210 nm