



1	Estriol 10 µg/mL	2	16α-Hydroxyestosterone 10 µg/mL	3	β ]Estradiol 10 µg/mL
4	Estrone 10 µg/mL	5	β ]Estradiol-3-methylester 10 µg/mL		

#### HPLC Conditions

Column	CAPCELL PAK C18 SG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	CH <sub>3</sub> CN / H <sub>2</sub> O = 55 / 45
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 210 nm