



1	Testosterone	2	$\beta$ ]Estradiol	3	4-Androstene-3,17-dione
	100 ppm		100 ppm		100 ppm
4	Estrone				
	100 ppm				

#### HPLC Conditions

Column	CAPCELL PAK CN UG120 S5 : 4.6 mm i.d.×250 mm
Mobile phase	H <sub>2</sub> O / CH <sub>3</sub> CN = 65 / 35
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 220 nm
Inj. vol.	5 $\mu$ L