



- | | | |
|---------------|-----------------|------------------|
| 1 DOPA | 2 Creatinine | 3 Norepinephrine |
| 4 Epinephrine | 5 Octopamine | 6 Normetanephine |
| 7 Dopamine | 8 Isoproterenol | 9 Metanephine |
| 10 Tyramine | 11 Serotonin | |

HPLC Conditions

Column CAPCELL PAK SCX UG80 : 4.6 mm i.d.×250 mm
Mobile phase 20mmol/L CH₃COONa, 0.02mol/L CH₃COOH (pH5.61)
Flow rate 1.0 mL/min
Temperature 40 °C
Detect UV 254 nm