



1	Teststerone	2	Androstenedione	3	17 α -Hydroxyprogesterone
	2.5 pmol		2.5 pmol		2.5 pmol
4	Progesterone				
	2.5 pmol				

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S3 : 2.0 mm i.d.×100 mm
Mobile phase	H ₂ O / CH ₃ CN = 68 / 32
Flow rate	300 μ L/min
Temperature	40 °C
Detect	UV 245 nm
Inj. vol.	2 μ L