



- | | | | | | |
|---|----------------------------------|---|---------------------------------|---|---------------------------------|
| 1 | 6 α -Hydroxytestosterone | 2 | 7 α -Hydroxytestosterone | 3 | 6 β -Hydroxytestosterone |
| 4 | 16 α -Hydroxytestosterone | 5 | 16 β -Hydroxytestosterone | 6 | 2 α -Hydroxytestosterone |
| 7 | 2 β -Hydroxytestosterone | 8 | Androsterone | 9 | Methyltestosterone |

HPLC Conditions

Column	CAPCELL PAK C18 UG120 S3 : 4.6 mm i.d.×150 mm
Mobile phase	H ₂ O / CH ₃ OH / THF = 55 / 38 / 7
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm