



1 Tryptic digest of casein

•^a%_ođ”½%_ož’O, ðfJf[fCf“”Z“x, Å5mg/mL

HPLC Conditions

Column CAPCELL PAK C18 ACR S5 : 1.0 mm i.d.×150 mm

Mobile phase (A) 0.05vol% TFA, H₂O
(B) 0.05vol% TFA, CH₃CN
B 10% (0min) - 50% (60min) Gradient

Flow rate 50 µL/min

Temperature 40 °C

Detect UV 210 nm

Inj. vol. 1 µL