



1 Tryptic digest of casein

•<sup>a</sup>%<sub>o</sub>đ”½%<sub>o</sub>ž‘O, ÌfJf[fCf“”Z“‘x, Å5mg/mL

#### HPLC Conditions

Column	CAPCELL PAK C18 UG120 S5 : 1.0 mm i.d.×150 mm
Mobile phase	(A) 0.05vol%TFA, H <sub>2</sub> O (B) 0.05vol%TFA, CH <sub>3</sub> CN B 10% (0min) - 50% (60min) Gradient
Flow rate	50 µL/min
Temperature	40 °C
Detect	UV 210 nm
Inj. vol.	1 µL