



1 Tryptic digest of myoglobin

•^a%_ođ”½%_ož’O, ìf~fIfOf frf“”Z“x, Å7.5mg/mL

HPLC Conditions

Column	CAPCELL PAK C18 ACR S5 : 1.0 mm i.d.×150 mm
Mobile phase	(A) 0.05vol% TFA, H ₂ O (B) 0.05vol% TFA, CH ₃ CN B 10% (0min) - 50% (60min) Gradient
Flow rate	50 µL/min
Temperature	40 °C
Detect	UV 210 nm
Inj. vol.	1 µL