



1	Thiuram	2	Iprodione	3	Bensulide
	0.2 µg/mL		0.2 µg/mL		0.5 µg/mL

#### HPLC Conditions

Column	CAPCELL PAK C <sub>18</sub> SG120 S5 : 4.6 mm i.d.×150 mm
Mobile phase	50mmol/L KH <sub>2</sub> PO <sub>4</sub> (pH3.3) / CH <sub>3</sub> CN = 50 / 50
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 230 nm
Inj. vol.	20 µL