



1 Oxine-copper  
0.02 ppm

2 Asulam  
0.1 ppm

#### HPLC Conditions

Column CAPCELL PAK C18 UG120 S5 : 1.0 mm i.d.×250 mm  
Mobile phase 2vol% H<sub>3</sub>PO<sub>4</sub> / CH<sub>3</sub>CN = 95 / 5  
Flow rate 50 μL/min  
Temperature 40 °C  
Detect UV 250 nm (0-8min), 270 nm (8-15min)  
Inj. vol. 8 μL