



1	<i>p</i> -Formylbenzoic acid	2	Benzoic acid	3	<i>p</i> -Toluic acid
	25 ppm		25 ppm		25 ppm
4	<i>p</i> -Etylbenzoic acid				
	25 ppm				

HPLC Conditions

Column	CAPCELL PAK C ₁₈ UG120 S5 : 1.5 mm i.d.×150 mm
Mobile phase	0.1vol% Phosphate, H ₂ O / CH ₃ CN = 75 / 25
Flow rate	100 μL/min
Temperature	R.T.
Detect	UV 247 nm
Inj. vol.	1 μL