



| | | | | | |
|---|---------------|---|--------------|---|-------------------------|
| 1 | Vanillic acid | 2 | Ferulic acid | 3 | <i>o</i> -Coumaric acid |
| | 100 ppm | | 100 ppm | | 100 ppm |

HPLC Conditions

| | |
|--------------|--|
| Column | CAPCELL PAK C ₁₈ UG120 S5 : 4.6 mm i.d.×250 mm |
| Mobile phase | 0.1mol/L H ₃ PO ₄ / CH ₃ CN = 75 / 25 |
| Flow rate | 1.0 mL/min |
| Temperature | R.T. |
| Detect | UV 254 nm |
| Inj. vol. | 5 μL |