



1	プロピオン酸 Propionic acid	2	酢酸 Acetic acid	3	コハク酸 Succinic acid
4	乳酸 Lactic acid	5	ピログルタミン酸 Pyroglutamic acid		

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	5mmol/L KH ₂ PO ₄ (pH4.0, H ₃ PO ₄), CH ₃ CN / H ₂ O = 70 / 30
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 210 nm