



1	レゾルシノール	2	スウェルチアマリン	3	<i>p</i> -ヒドロキシ安息香酸
	Resorcinol		Swertiamarin		<i>p</i> -Hydroxybenzoic acid
	1000 ppm		200 ppm		100 ppm

HPLC Conditions	
Column	CAPCELL PAK C ₁₈ MG S5 : 4.6 mm i.d. × 250 mm
Mobile phase	H ₂ O / CH ₃ CN = 87 / 13 (v/v), pH2.5 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 238 nm
Inj. vol.	5 μL