



1	L-リンゴ酸	2	マロン酸	3	乳酸
	L-Malic acid		Malonic acid		Lactic acid
	500 ppm		500 ppm		500 ppm
4	酢酸	5	クエン酸	6	コハク酸
	Acetic acid		Citric acid		Succinic acid
	500 ppm		500 ppm		500 ppm
7	プロピオン酸				
	Propionic acid				
	500 ppm				

HPLC Conditions

Column	CAPCELL PAK C ₁₈ MG S5 : 4.6 mm i.d. × 150 mm
Mobile phase	0.1vol% H ₃ PO ₄ , H ₂ O / CH ₃ CN = 97.5 / 2.5 (v/v)
Flow rate	1.0 mL/min
Temperature	35
Detect	UV 210 nm
Inj. vol.	7 μL

