



1	グリコール酸 Glycolic acid	2	ギ酸 Formic acid	3	酢酸 Acetic acid
---	-------------------------	---	-------------------	---	-------------------

HPLC Conditions	
Column	CAPCELL PAK C <sub>18</sub> MG S5 : 2.0 mm i.d. × 250 mm
Mobile phase	0.1vol% H <sub>3</sub> PO <sub>4</sub> , H <sub>2</sub> O / CH <sub>3</sub> CN = 97.5 / 2.5 (v/v)
Flow rate	200 μL/min
Temperature	35
Detect	UV 210 nm