



1	ビオチン (VH)	2	リボフラビン (VB ₂)	3	ニコチン酸
	Biotin (VH)		Riboflavin (VB ₂)		Nicotinic acid
4	塩酸ピリドキシン (VB ₆)	5	葉酸 (VBc)	6	ニコチンアミド (PP)
	Pyridoxine hydrochloride (VB ₆)		Folic acid (VBc)		Nicotinamide (PP)

HPLC Conditions

Column	CAPCELL PAK NH ₂ UG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	CH ₃ CN / 25mmol/L KH ₂ PO ₄ = 80 / 20, pH2.0 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 220 nm

Inj. vol.

20 μ L