



1	ビオチン (VH)	2	リボフラビン (VB ₂)	3	アスコルビン酸 (VC)
	Biotin (VH)		Riboflavin (VB ₂)		Ascorbic acid (VC)
	80 µg/mL		80 µg/mL		100 µg/mL
4	ニコチン酸	5	葉酸 (VBc)	6	ニコチンアミド (PP)
	Nicotinic acid		Folic acid (VBc)		Nicotinamide (PP)
	20 µg/mL		160 µg/mL		20 µg/mL

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	CH ₃ CN / 25mmol/L KH ₂ PO ₄ = 80 / 20, pH2.0 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 210 nm

Inj. vol.

10 μ L