



| | | | |
|---|------------------|---|--------------------|
| 1 | イソアスコルビン酸 | 2 | アスコルビン酸 (VC) |
| | Isoascorbic acid | | Ascorbic acid (VC) |

HPLC Conditions

| | |
|--------------|--|
| Column | CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d. × 250 mm |
| Mobile phase | 2mmol/L KH ₂ PO ₄ , CH ₃ CN / H ₂ O = 85 / 15, pH2.0 (H ₃ PO ₄) |
| Flow rate | 1.0 mL/min |
| Temperature | 40 |
| Detect | UV 254 nm |