



1	ビタミンK ₃ Vitamin K ₃ 24 ppm	2	ビタミンA Vitamin A 320 ppm	3	ビタミンD ₂ Vitamin D ₂ 500 ppm
4	ビタミンD ₃ Vitamin D ₃ 100 ppm	5	ビタミンE Vitamin E 500 ppm	6	ビタミンEアセテート Vitamin E acetate 5000 ppm
7	ビタミンK ₁ Vitamin K ₁ 200 ppm				

HPLC Conditions

Column	CAPCELL PAK C ₁₈ MG S5 : 4.6 mm i.d. × 250 mm
Mobile phase	CH ₃ CN / CH ₃ OH = 80 / 20
Flow rate	1.0 mL/min
Temperature	35

Detect	UV 254 nm
Inj. vol.	10 μ L