



1	コンドロΔDi-0S	2	コンドロΔDi-6S	3	コンドロΔDi-4S
	Chondro ΔDi-0S		Chondro ΔDi-6S		Chondro ΔDi-4S
	0.6 μg/mL		0.6 μg/mL		0.6 μg/mL

HPLC Conditions

Column	CAPCELL PAK NH ₂ UG80 S5 : 1.5 mm i.d. × 250 mm
Mobile phase	0.25mol/L KH ₂ PO ₄
Flow rate	100 μL/min
Temperature	40
Detect	UV 232 nm
Inj. vol.	10 μL