



1	5'-AMP	2	5'-ADP	3	5'-ATP
	5'-AMP		5'-ADP		5'-ATP
	100 µg/mL		100 µg/mL		100 µg/mL

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	1.0mol/L (NH ₄) ₂ HPO ₄ (pH2.0)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 260 nm
Inj. vol.	10 µL