



1	5'-CMP	2	5'-AMP	3	5'-UMP
	5'-CMP		5'-AMP		5'-UMP
	100 ppm		100 ppm		100 ppm
4	5'-IMP	5	5'-GMP		
	5'-IMP		5'-GMP		
	100 ppm		100 ppm		

## HPLC Conditions

Column	CAPCELL PAK NH <sub>2</sub> UG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	50mmol/L (NH <sub>4</sub> ) <sub>2</sub> HPO <sub>4</sub> (pH3.0)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 254 nm
Inj. vol.	10 μL

