



1	5'-CMP	2	5'-AMP	3	5'-UMP
	5'-CMP		5'-AMP		5'-UMP
	100 µg/mL		100 µg/mL		100 µg/mL
4	5'-IMP	5	5'-GMP		
	5'-IMP		5'-GMP		
	100 µg/mL		100 µg/mL		

HPLC Conditions

Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	50mmol/L (NH ₄) ₂ HPO ₄ (pH3.0)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 260 nm
Inj. vol.	10 µL

