



1	ウラシル	2	チミン	3	グアニン
	Uracil		Thymine		Guanine
	5 µg/mL		5 µg/mL		15 µg/mL
4	シトシン	5	アデニン		
	Cytosine		Adenine		
	30 µg/mL		15 µg/mL		

## HPLC Conditions

Column	CAPCELL PAK SCX UG80 : 4.6 mm i.d. × 250 mm
Mobile phase	200mmol/L NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub> (pH3.50)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 254 nm
Inj. vol.	10 µL