



1	ホルムアルデヒド-2,4-DNPH	2	アセトアルデヒド-2,4-DNPH	3	プロピオンアルデヒド-2,4-DNPH
	Formaldehyde-2,4-DNPH		Acetaldehyde-2,4-DNPH		Propionaldehyde-2,4-DNPH
	1 ppm		1 ppm		1 ppm
4	イソブチルアルデヒド-2,4-DNPH	5	吉草酸アルデヒド-2,4-DNPH	6	イソ吉草酸アルデヒド-2,4-DNPH
	Isobutyraldehyde-2,4-DNPH		Valeraldehyde-2,4-DNPH		isovaleraldehyde-2,4-DNPH
	1 ppm		1 ppm		1 ppm

HPLC Conditions

Column	CAPCELL PAK C ₁₈ MG S5 : 1.5 mm i.d. × 150 mm
Mobile phase	CH ₃ CN / H ₂ O = 60 / 40
Flow rate	100 μL/min
Temperature	40
Detect	UV 360 nm
Inj. vol.	0.5 μL