



1	2-プロピルピリジン 2-Propylpyridine	2	2-エチルピリジン 2-Ethylpyridine	3	2-メチルピリジン 2-Methylpyridine
4	ピリジン Pyridine				

HPLC Conditions	
Column	CAPCELL PAK NH ₂ SG80 S5 : 4.6 mm i.d. × 250 mm
Mobile phase	5mmol/L KH ₂ PO ₄ , CH ₃ CN / H ₂ O = 80 / 20, pH2.0 (H ₃ PO ₄)
Flow rate	1.0 mL/min
Temperature	40
Detect	UV 254 nm