



1	<i>p</i> -クレゾール	2	<i>m</i> -クレゾール	3	<i>o</i> -クレゾール
	<i>p</i> -cresol		<i>m</i> -cresol		<i>o</i> -cresol
	50 µg/mL		50 µg/mL		50 µg/mL

HPLC Conditions

Column	CAPCELL PAK C ₁₈ MG S3 : 1.5 mm i.d. × 75 mm
Mobile phase	10mmol/L β-Cyclodextrin, 50mmol/L Sodium phosphate buffer (pH5.0) / CH ₃ CN = 90 / 10
Flow rate	200 µL/min
Temperature	40
Detect	UV 277 nm
Inj. vol.	2.0 µL