



1	フマル酸 Fumaric acid	2	アクリル酸 Acrylic acid	3	イタコン酸 Itaconic acid
4	<i>p</i> -ヒドロキシ安息香酸 <i>p</i> -Hydroxybenzoic acid	5	バニリン酸 Vanillic acid		

HPLC Conditions

Column CAPCELL PAK C₁₈ UG120 S5 : 4.6 mm i.d. × 250 mm
 CAPCELL PAK C₁₈ MG S5 : 4.6 mm i.d. × 250 mm
 CAPCELL PAK C₁₈ AQ S5 : 4.6 mm i.d. × 250 mm

Mobile phase 0.1vol% H₃PO₄ / CH₃OH = 90 / 10

Detect UV 210 nm