

製品ラインナップ 追加のご案内

< 3 μm 高耐圧化ラインナップの追加 >

粒子径3 μmで50MPa耐圧仕様のカラムをラインナップに追加しました。

製品番号	製品名	粒子径 (μm)	内径 (mm)	長さ (mm)	標準価格 (税抜)
92499	CAPCELL PAK C18 MGII-H (S3) 4.6X250	3	4.6	250	¥65,000
92498	CAPCELL PAK C18 MGII-H (S3) 4.6X150	3	4.6	150	¥60,000
92497	CAPCELL PAK C18 MGII-H (S3) 4.6X100	3	4.6	100	¥52,000
92432	CAPCELL PAK C18 MGII-H (S3) 2.0X150	3	2.0	150	¥76,000
92431	CAPCELL PAK C18 MGII-H (S3) 2.0X100	3	2.0	100	¥71,000
92430	CAPCELL PAK C18 MGII-H (S3) 4.0X50	3	2.0	50	¥66,000
90890	CAPCELL PAK C18 MG-H (S3) 4.6X250	3	4.6	250	¥62,000
90889	CAPCELL PAK C18 MG-H (S3) 4.6X150	3	4.6	150	¥58,000
90888	CAPCELL PAK C18 MG-H (S3) 4.6X100	3	4.6	100	¥50,000
90825	CAPCELL PAK C18 MG-H (S3) 2.0X150	3	2.0	150	¥73,000
90824	CAPCELL PAK C18 MG-H (S3) 2.0X100	3	2.0	100	¥68,000
90823	CAPCELL PAK C18 MG-H (S3) 2.0X50	3	2.0	50	¥63,000
93349	CAPCELL PAK ADME-HR-H (S3) 4.6X250	3	4.6	250	¥96,000
93329	CAPCELL PAK ADME-HR-H (S3) 2.1X250	3	2.1	250	101,000
93328	CAPCELL PAK ADME-HR-H (S3) 2.1X150	3	2.1	150	85,000
93327	CAPCELL PAK ADME-HR-H (S3) 2.1X100	3	2.1	100	78,000
93319	CAPCELL PAK ADME-HR-H (S3) 2.1X50	3	2.1	50	79,000
92399	CAPCELL PAK C8 DD-H (S3) 4.6X250	3	4.6	250	¥65,000
92342	CAPCELL PAK C8 DD-H (S3) 2.0X150	3	2.0	150	¥68,000
92341	CAPCELL PAK C8 DD-H (S3) 2.0X100	3	2.0	100	¥68,000
92340	CAPCELL PAK C8 DD-H (S3) 2.0X50	3	2.0	50	¥63,000
92170	CAPCELL PAK C18 AQ-H (S3) 4.6X250	3	4.6	250	¥78,000
92139	CAPCELL PAK C18 AQ-H (S-3) 2.0X150	3	2.0	150	¥70,000
92138	CAPCELL PAK C18 AQ-H (S-3) 2.0X100	3	2.0	100	¥68,000
92137	CAPCELL PAK C18 AQ-H (S-3) 2.0X50	3	2.0	50	¥63,000
92260	CAPCELL PAK C18 ACR-H (S3) 4.6X250	3	4.6	250	¥65,000
92239	CAPCELL PAK C18 ACR-H (S-3) 2.0X150	3	2.0	150	¥68,000
92238	CAPCELL PAK C18 ACR-H (S-3) 2.0X100	3	2.0	100	¥68,000
92237	CAPCELL PAK C18 ACR-H (S-3) 2.0X50	3	2.0	50	¥63,000
44560	SUPERIOREX II ODS-H (S3) 4.6X250	3	4.6	250	¥69,000
90889	CAPCELL PAK C18 MG-H (S-3) 4.6X150	3	4.6	150	¥58,000
90888	CAPCELL PAK C18 MG-H (S-3) 4.6X100	3	4.6	100	¥50,000
90825	CAPCELL PAK C18 MG-H (S-3) 2.0X150	3	2.0	150	¥73,000
90824	CAPCELL PAK C18 MG-H (S-3) 2.0X100	3	2.0	100	¥68,000
90823	CAPCELL PAK C18 MG-H (S-3) 2.0X50	3	2.0	50	¥63,000
61308	CAPCELL PAK C18 UG120-H (S-3) 2.0X150	3	2.0	150	¥73,000
61306	CAPCELL PAK C18 UG120-H (S-3) 2.0X100	3	2.0	100	¥68,000
61307	CAPCELL PAK C18 UG120-H (S-3) 2.0X50	3	2.0	50	¥63,000

< CAPCELL PAK SEC > サイズ排除モードカラム

グレード	細孔径 (Å)	粒子径 (μm)	分画範囲	耐圧 (MPa)	使用pH範囲
SEC 125	125	2/5/10	5 – 150k Da (タンパク質)	35/8/2	2-7.5
SEC 250	250	2/5/10	10 – 500k Da (タンパク質)	35/8/2	2-7.5
SEC 1000	1000	5	100 – 7,500k Da (タンパク質), 30 – 3,000 nt (ssRNA)	8	2-7.5

製品番号	製品名	粒子径 (μm)	内径 (mm)	長さ (mm)	標準価格 (税抜)
82111	CAPCELL PAK SEC125 (S2) 4.6X150	2	4.6	150	¥190,000
82112	CAPCELL PAK SEC125 (S2) 4.6X300	2	4.6	300	¥270,000
82135	CAPCELL PAK SEC125 (S5) 7.8X300	5	7.8	300	¥180,000
82145	CAPCELL PAK SEC125 (S10) 7.8X300	10	7.8	300	¥120,000
82211	CAPCELL PAK SEC250 (S2) 4.6X150	2	4.6	150	¥190,000
82212	CAPCELL PAK SEC250 (S2) 4.6X300	2	4.6	300	¥270,000
82235	CAPCELL PAK SEC250 (S5) 7.8X300	5	7.8	300	¥180,000
82245	CAPCELL PAK SEC250 (S10) 7.8X300	10	7.8	300	¥120,000
82535	CAPCELL PAK SEC1000 (S5) 7.8X300	5	7.8	300	¥230,000